

# HUNK Polyclonal Antibody

Cat No: HR1AP6842

For research use only

## Overview

Product Name	HUNK Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	HUNK
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HUNK; MAKV; Hormonally up-regulated neu tumor-associated kinase; B19; Serine/threonine-protein kinase MAK-V
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. SNF1 subfamily.,similarity:Contains 1 protein kinase domain.,