

DGK-? Polyclonal Antibody

Cat No: HR1AP5925

For research use only

Overview

Product Name	DGK-? Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	142kDa
GeneID?Human?	DGKK
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DGKK; Diacylglycerol kinase kappa; DAG kinase kappa; 142 kDa diacylglycerol kinase; Diglyceride kinase kappa; DGK-kappa
Background	diacylglycerol kinase kappa(DGKK) Homo sapiens The protein encoded by this gene is an enzyme that phosphorylates diacylglycerol, converting it to phosphatidic acid. The encoded protein is a membrane protein and is inhibited by hydrogen peroxide. Variations in this gene have been associated with hypospadias. [provided by RefSeq, Mar 2011],