

## Dyrk4 Polyclonal Antibody

Cat No: HR1AP8990

For research use only

### Overview

Product Name	Dyrk4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	DYRK4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DYRK4; Dual specificity tyrosine-phosphorylation-regulated kinase 4
Background	<p>dual specificity tyrosine phosphorylation regulated kinase 4(DYRK4) Homo sapiens This gene encodes an enzyme that belongs to a conserved family of serine/threonine protein kinases. Members of this dual specificity kinase family are thought to function in the regulation of cell differentiation and proliferation, survival, and in development. Alternate splicing results in multiple transcript variants. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Aug 2013],</p>