

Alpha-protein Kinase 2 Polyclonal Antibody

Cat No: HR1AP5768

For research use only

Overview

Product Name	Alpha-protein Kinase 2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	ALPK2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ALPK2; HAK; Alpha-protein kinase 2; Heart alpha-protein kinase
Background	function:Kinases that recognize phosphorylation sites in which the surrounding peptides have an alpha-helical conformation.,similarity:Belongs to the protein kinase superfamily. Alpha-type protein kinase family. ALPK subfamily.,similarity:Contains 1 alpha-type protein kinase domain.,similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,