

ADCK2 Polyclonal Antibody

Cat No: HR1AP9068

For research use only

Overview

Product Name	ADCK2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	69kDa
GeneID?Human?	ADCK2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ADCK2; AARF; Uncharacterized aarF domain-containing protein kinase 2
Background	function:The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).,sequence caution:Aberrant splicing.,similarity:Belongs to the protein kinase superfamily. ADCK protein kinase family.,similarity:Contains 1 protein kinase domain.,