

ACYP2 Polyclonal Antibody

Cat No: HR1AP9211

For research use only

Overview

Product Name	ACYP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	ACYP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ACYP2; ACYP; Acylphosphatase-2; Acylphosphatase; muscle type isozyme; Acylphosphate phosphohydrolase 2
Background	<p>acylphosphatase 2(ACYP2) Homo sapiens Acylphosphatase can hydrolyze the phosphoenzyme intermediate of different membrane pumps, particularly the Ca²⁺/Mg²⁺-ATPase from sarcoplasmic reticulum of skeletal muscle. Two isoenzymes have been isolated, called muscle acylphosphatase and erythrocyte acylphosphatase on the basis of their tissue localization. This gene encodes the muscle-type isoform (MT). An increase of the MT isoform is associated with muscle differentiation. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2016].</p>