

PDP2 Polyclonal Antibody

Cat No: HR1AP4915

For research use only

Overview

Product Name	PDP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	59kDa
GeneID?Human?	PDP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PDP2; KIAA1348; [Pyruvate dehydrogenase [acetyl-transferring]]-phosphatase 2, mitochondrial; PDP 2; Pyruvate dehydrogenase phosphatase catalytic subunit 2; PDPC 2
Background	pyruvate dehydrogenase phosphatase catalytic subunit 2(PDP2) Homo sapiens This gene is a mitochondrial protein that functions as a phosphatase and is involved in the enzymatic resetting of the pyruvate dehydrogenase complex. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Aug 2016],