

## ?-protein Kinase 1 Polyclonal Antibody

Cat No: HR1AP4857

For research use only

### Overview

Product Name	?-protein Kinase 1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	139kDa
GeneID?Human?	ALPK1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ALPK1; KIAA1527; LAK; Alpha-protein kinase 1; Chromosome 4 kinase; Lymphocyte alpha-protein kinase
Background	alpha kinase 1(ALPK1) Homo sapiens This gene encodes an alpha kinase. Mice which were homozygous for disrupted copies of this gene exhibited coordination defects (PMID: 21208416). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011],