

ULK3 Polyclonal Antibody

Cat No: HR1AP4785

For research use only

Overview

Product Name	ULK3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	ULK3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ULK3; Serine/threonine-protein kinase ULK3; Unc-51-like kinase 3
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family. APG1/unc-51/ULK1 subfamily.,similarity:Contains 1 protein kinase domain.,