

PPP4R1L Polyclonal Antibody

Cat No: HR1AP4365

For research use only

Overview

Product Name	PPP4R1L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	PPP4R1L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PPP4R1L; C20orf192; Putative serine/threonine-protein phosphatase 4 regulatory subunit 1-like
Background	function:May be a regulatory subunit of serine/threonine-protein phosphatase 4.,similarity:Contains 4 HEAT repeats.,