

## PDXDC1 Polyclonal Antibody

Cat No: HR1AP4283

For research use only

### Overview

Product Name	PDXDC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	87kDa
GeneID?Human?	PDXDC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PDXDC1; KIAA0251; Pyridoxal-dependent decarboxylase domain-containing protein 1
Background	cofactor:Pyridoxal phosphate.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the group II decarboxylase family.,