

GAP1-InsP4 BP Polyclonal Antibody

Cat No: HR1AP3508

For research use only

Overview

Product Name	GAP1-InsP4 BP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	RASA3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RASA3; Ras GTPase-activating protein 3; GAP1(IP4BP); Ins P4-binding protein
Background	RAS p21 protein activator 3(RASA3) Homo sapiens This gene encodes a protein that binds inositol 1,3,4,5-tetrakisphosphate and stimulates the GTPase activity of Ras p21. This protein functions as a negative regulator of the Ras signalling pathway. It is localized to the cell membrane via a pleckstrin homology (PH) domain in the C-terminal region. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016],