

ADCK5 Polyclonal Antibody

Cat No: HR1AP2714

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | ADCK5 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IF,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 66kDa |
| GeneID?Human? | ADCK5 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | ADCK5; Uncharacterized aarF domain-containing protein kinase 5 |
| Background | function:The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).,similarity:Belongs to the protein kinase superfamily. ADCK protein kinase family.,similarity:Contains 1 protein kinase domain., |