

IRAK-2 Polyclonal Antibody

Cat No: HR1AP3753

For research use only

Overview

Product Name	IRAK-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	IRAK2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IRAK2; Interleukin-1 receptor-associated kinase-like 2; IRAK-2
Background	interleukin 1 receptor associated kinase 2(IRAK2) Homo sapiens IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB. [provided by RefSeq, Jul 2008],