

JAK1 Rabbit Polyclonal Antibody

Cat No: HR1AP9673

For research use only

Overview

Product Name	JAK1 Rabbit Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	130kDa
GeneID?Human?	JAK1 JAK1A JAK1B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Tyrosine-protein kinase JAK1 (EC 2.7.10.2) (Janus kinase 1) (JAK-1)
Background	Janus kinase 1(JAK1) Homo sapiens This gene encodes a membrane protein that is a member of a class of protein-tyrosine kinases (PTK) characterized by the presence of a second phosphotransferase-related domain immediately N-terminal to the PTK domain. The encoded kinase phosphorylates STAT proteins (signal transducers and activators of transcription) and plays a key role in interferon-alpha/beta and interferon-gamma signal transduction. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016].