

HXK III Polyclonal Antibody

Cat No: HR1AP6851

For research use only

Overview

Product Name	HXK III Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	99kDa
GeneID?Human?	HK3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HK3; Hexokinase-3; Hexokinase type III; HK III
Background	hexokinase 3(HK3) Homo sapiens Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. This gene encodes hexokinase 3. Similar to hexokinases 1 and 2, this allosteric enzyme is inhibited by its product glucose-6-phosphate. [provided by RefSeq, Apr 2009],