

DAPK1 (phospho Ser308) Polyclonal Antibody

Cat No: HR1AP6075

For research use only

Overview

Product Name	DAPK1 (phospho Ser308) Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	DAPK1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DAPK1; DAPK; Death-associated protein kinase 1; DAP kinase 1
Background	death associated protein kinase 1(DAPK1) Homo sapiens Death-associated protein kinase 1 is a positive mediator of gamma-interferon induced programmed cell death. DAPK1 encodes a structurally unique 160-kD calmodulin dependent serine-threonine kinase that carries 8 ankyrin repeats and 2 putative P-loop consensus sites. It is a tumor suppressor candidate. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013],