

Blk (phospho Tyr501) Polyclonal Antibody

Cat No: HR1AP8255

For research use only

Overview

Product Name	Blk (phospho Tyr501) Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	BLK
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BLK; Tyrosine-protein kinase Blk; B lymphocyte kinase; p55-Blk
Background	<p>BLK proto-oncogene, Src family tyrosine kinase(BLK) Homo sapiens This gene encodes a nonreceptor tyrosine-kinase of the src family of proto-oncogenes that are typically involved in cell proliferation and differentiation. The protein has a role in B-cell receptor signaling and B-cell development. The protein also stimulates insulin synthesis and secretion in response to glucose and enhances the expression of several pancreatic beta-cell transcription factors. [provided by RefSeq, Aug 2010].</p>