

# IKK $\gamma$ (phospho Ser85) Polyclonal Antibody

Cat No: HR1AP8912

For research use only

## Overview

Product Name	IKK $\gamma$ (phospho Ser85) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	IKBKG
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IKBKG; FIP3; NEMO; NF-kappa-B essential modulator; NEMO; FIP-3; I $\kappa$ B kinase-associated protein 1; IKKAP1; Inhibitor of nuclear factor kappa-B kinase subunit gamma; I-kappa-B kinase subunit gamma; IKK-g
Background	inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase gamma(IKBKG) Homo sapiens This gene encodes the regulatory subunit of the inhibitor of kappaB kinase (IKK) complex, which activates NF-kappaB resulting in activation of genes involved in inflammation, immunity, cell survival, and other pathways. Mutations in this gene result in incontinentia pigmenti, hypohidrotic ectodermal dysplasia, and several other types of immunodeficiencies. A pseudogene highly similar to this locus is located in an adjacent region of the X chromosome. [provided by RefSeq, Mar 2016],