

## PP4R2 Polyclonal Antibody

Cat No: HR1AP6023

For research use only

### Overview

Product Name	PP4R2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	PPP4R2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PPP4R2; SBBI57; Serine/threonine-protein phosphatase 4 regulatory subunit 2
Background	protein phosphatase 4 regulatory subunit 2(PPP4R2) Homo sapiens The protein encoded by this gene is a regulatory subunit of the serine/threonine-protein phosphatase 4 complex. In addition to being required for efficient DNA double strand break repair, this complex plays a role in organization of microtubules at centrosomes and processing of spliceosomal snRNPs. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2015],