



# Kit(Density Reagent)stein Cholesterol (LDL-C) Colorimetric Assay

Cat No: HR3BC1148

For research use only

## Overview

Detection Method	Colorimetric method
Storage	2-8?
Instrument	Spectrophotometer(546 nm)
Assay Time	
Validity	12
Assay Type	Quantitative
Sample Type	Serum, plasma, cells, culture supernatant, tissue
Synonyms	LDL-C
Instrument	Spectrophotometer(546 nm)
Detection Principle	The coloured substance have a maximum absorption peak at 546 nm. Measure the OD value at 546 nm and the LDL-C content in the sample can be calculated.
Reagents	
Labware	
Size	100Assays
Sensitivity	
Detection Range	0.2-12 mmol/L
Recovery Rate	