

Recombinant Human NIK (N-GST tag)

Cat No:HR2R1792

For research use only

Overview

Quantity	10 ?g
Gene Symbol	NIK
Gene ID	N/A
Accession	NM_003954
Alternative Name	MAP3K14, HS, HSNIK, FTDCR1B
Species	Human
Source	Insect cells
Description	<p>NIK is a mitogen-activated protein kinase kinase kinase 14 (MAP3K14), which binds to TRAF2 and stimulates NF-kappaB activity. NIK shares sequence similarity with several other MAPKK kinases and participates in NF-kappaB-inducing signalling cascade common to receptors of the tumour-necrosis/nerve-growth factor (TNF/NGF) family and to the interleukin-1 type-I receptor. NIK is expressed in primary human cells and in inflamed rheumatoid arthritis tissue and plays a selective role in signaling by the lymphotoxin-beta receptor . NIK is a therapeutic target in the immune and bone-destructive components of inflammatory arthritis.</p>
Functions	The specific activity of NIK was determined to be 8 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	108
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.